

IMPORTANT INFORMATION TO READ and RETURN

Installation Requirements for a Whitley A155 HEPA Anaerobic Workstation

Thank you for choosing one of our products for your laboratory. To enable our engineers to perform an efficient, trouble-free installation please study, complete and return this form to us by fax: **01274 531197** or email: **service@dwscientific.co.uk**. Should you have any questions, please do not hesitate to contact us, as we are here to help. When we have received the completed form, our Service Department will contact you to arrange a mutually convenient installation date.

The following information represents the ideal requirement.
Please contact us IMMEDIATELY if your intended location does not match this specification.



Access Requirements

Please note that this product is **wider than a standard door** so please ensure that you consider access carefully.

☐

Space Requirements

The weight of the equipment is 350 kg. This workstation comes complete with a bespoke trolley. Please ensure that sufficient floor space is allocated. If you wish to locate the workstation on a bench, please contact us to discuss the implications.

External Dimensions

| Width | Depth | Height | Height Including Trolley |
|-------|-------|--------|--------------------------|
| mm | mm | mm | mm |
| 2875 | 1056 | 993 | 1911 |

☐

Please Note:

If the **Heat Removal System** is fitted, please allow a clearance of 600mm behind the workstation.
Some of this clearance is necessary for air circulation.

Gas Requirements

The incoming gas supplies must be terminated near the right hand side of the main chamber and fitted with leak-proof taps and pressure gauges.

The gas lines to which the equipment is attached are the responsibility of the user and should be constructed, tested and maintained to the standards specified within the British Compressed Gases Association (BCGA) Code of Practice CP4 (or international equivalent). Gas lines previously used for flammable gases must be purged prior to re-use.

Regulators should be fitted in accordance with the information contained in the table below and the various pressures strictly adhered to.

| Gas Type | Connection Details | Cylinder Regulator Required | Regulator Outlet Range | Flow Rate |
|-------------------|---|-----------------------------|------------------------|---|
| Anaerobic Gas Mix | ¼" BSP male fitting or connection for 6mm Polyurethane tubing. Hydrogen/Anaerobic Gas Mixture Regulator – Two Stage – order Code A01745 | Two stage | 4 - 6 bar (60-90 psi) | 40 litres per minute (dynamic) |
| Nitrogen | ¼" BSP male fitting or connection for 8mm Polyurethane tubing Nitrogen Regulator – Two Stage – order Code A01748 | Two Stage | 4 - 6 bar (60-90 psi) | Minimum 150 litres per minute (dynamic) |

Note: If the instrument is to be connected to a single cylinder of anaerobic gas mixture only, then this must be capable of delivering 150 litres/minute (sleeved ports) or 100 litres/minute (Instant Access Port).


Connection Type



Push in connection

☐

OR



Push on Connection

☐

Mains Requirements

Electricity Supply

240 volts

Wall Socket

1 x Three Pin, 13 Amp. Minimum 6A rating

☐

Other Considerations

Although the workstation should be located in a well ventilated area, avoid close proximity to air conditioning system and draughts caused by windows and doors.

Remember, if you do not have the required regulators you can order them from Don Whitley Scientific::

- Hydrogen/Anaerobic Gas Mixture Regulator – Two Stage – order Code **A01745**
- Nitrogen Regulator – Two Stage – order Code **A01748**

☐

Notes

If an existing unit is being taken in part exchange or is being removed from the laboratory, it must be de-contaminated before DWS staff handle the unit. A certificate or signed letter confirming the unit has been decontaminated must be given to our engineer.

☐

In the UK, delivery and installation are free of charge (unless otherwise agreed). If our engineers are unable to install the unit and a return journey is necessary, **a charge may be made**. Export customers, please refer to your local distributor.

It is essential that this form is completed and returned, to avoid delay to your installation.

THANK YOU FOR THINKING WHITLEY

Signature

Title

Print Name

Establishment