

IMPORTANT INFORMATION TO READ and RETURN

Installation Requirements for the APS One Pourer Stacker

Thank you for choosing one of our products for your laboratory. To enable our engineers to perform an efficient, trouble-free installation please study, complete and fax this form to us on **01274 531197**. Should you have any questions, please do not hesitate to contact us, as we are here to help. When we have received the completed form, our Service Department will contact you to arrange a mutually convenient installation date.

**The following information represents the ideal requirement.
Please contact us IMMEDIATELY if your intended location does not match this specification.**



Access Requirements		<input type="checkbox"/>								
For access, the dimensions stated below should be taken into account when checking the size of lifts, stairwells, doorways, etc.										
Space Requirements										
The weight of the equipment is 60kg. The bench allocated must be flat, level and of sufficient size to support the base fully.										
The unit should be sited as close as possible (preferably to the left) of the Media Preparator.										
External Dimensions										
<table border="1" style="width: 100%; border-collapse: collapse; text-align: center;"> <thead> <tr> <th style="padding: 2px;">Depth</th> <th style="padding: 2px;">Length</th> <th style="padding: 2px;">Height</th> </tr> </thead> <tbody> <tr> <td style="padding: 2px;"><i>mm</i></td> <td style="padding: 2px;"><i>mm</i></td> <td style="padding: 2px;"><i>mm</i></td> </tr> <tr> <td style="padding: 2px;">640</td> <td style="padding: 2px;">720</td> <td style="padding: 2px;">900</td> </tr> </tbody> </table>	Depth	Length	Height	<i>mm</i>	<i>mm</i>	<i>mm</i>	640	720	900	<input type="checkbox"/>
Depth	Length	Height								
<i>mm</i>	<i>mm</i>	<i>mm</i>								
640	720	900								
Mains Requirements										
Electricity Supply	Socket	<input type="checkbox"/>								
230 volts +/-10%, 50-60Hz or 110 volts +/-10%, 50-60Hz	1 x Three pin, 13 amp									

Delivery and installation are free of charge (unless otherwise agreed). If our engineers are unable to install the unit and a return journey is necessary, **a charge may be made.**

It is essential that this form is completed and returned to avoid delay to your installation.

THANK YOU FOR THINKING WHITLEY

Signature	Title
<input style="width: 90%; height: 20px;" type="text"/>	<input style="width: 90%; height: 20px;" type="text"/>
Print Name	Establishment
<input style="width: 90%; height: 20px;" type="text"/>	<input style="width: 90%; height: 20px;" type="text"/>