

Whitley A45 Anaerobic Workstation - Removable Front/Porthole Options



Instant Access Porthole System

The entry level version of the Whitley A45 Anaerobic Workstation is equipped with 3 x Instant Access Porthole Systems.

This version is ideal for users wishing to enter the workstation frequently for brief working periods.

230V: Order A07500 / 115V: Order A07501



Instant Access Porthole System / Conventional Porthole System /

This version of the Whitley A45 Anaerobic Workstation is equipped with 1 x Instant Access Porthole System and 2 x Conventional Porthole Systems with sleeves and cuffs.

This version is ideal for users wishing to enter the workstation less frequently for longer working periods.

230V: Order A07540 / 115V: Order A07541



Instant Access Porthole System / Removable Front

This version of the Whitley A45 Anaerobic Workstation is equipped with 3 x Instant Access Porthole Systems and a removable front.

The removable front covers an opening which measures 800mm (31.5") wide x 398mm (15.5") high.

This version is ideal for users wishing to:

- regularly "deep clean" the working and incubation chamber;
- introduce small items of equipment
- add/remove incubation samples in bulk.
- enter the workstation frequently for brief working periods.

230V: Order A07530 / 115V: Order A07531



Instant Access Porthole System / Conventional Porthole System / Removable Front

This version of the Whitley A45 Anaerobic Workstation is equipped with 1 x Instant Access Porthole System, 2 x Conventional Porthole Systems with sleeves and cuffs and a removable front.

The removable front covers an opening which measures 800mm (31.5") wide x 398mm (15.5") high.

This version is ideal for users wishing to:

- regularly "deep clean" the working and incubation chamber;
- introduce small items of equipment
- add/remove incubation samples in bulk.
- enter the workstation less frequently for longer working periods.

230V: Order A07520 / 115V: Order A07521

Option 3 can be converted to Option 4 at a customer's premises at a future date, but additional components must be included during construction. For a 230V version order A07510, for a 115V version order A07511.

If/when conversion is required - order A07209 - a removable front fitted with 2 x Conventional Porthole Systems.